

Safety Data Sheet

according to WHS Regulations

Printing date 28.07.2023

Revision: 28.07.2023

1 Identification

- **Product identifier**
- **Trade name:** *DermaDam™ Synthetic*
- **Article number:** *SDS 163-001.04R02, 299, 330*
- **CAS Number:**
9003-31-0
- **Relevant identified uses of the substance or mixture and uses advised against** *Professional Dental Dam*
- **Application of the substance / the mixture** *Professional Dental Dam*
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
*Ultradent Products, Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
Ultradent Australia Pty Ltd.
Level 22/2 Market Street
Sydney NSW 2000
Australia
Email: info.anz@ultradent.com
Toll Free: 1-800-290929*
- **Further information obtainable from:** *Customer Service*
- **Emergency telephone number:**
*CHEMTREC (NORTH AMERICA) : (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887*

2 Hazard(s) Identification

- **Classification of the substance or mixture**
The substance is not classified, according to the Globally Harmonized System (GHS).
- **Additional information:**
*Health hazard: Polyisoprene is generally considered non-hazardous and is safe for use by persons who know they are susceptible to natural latex allergies. However, the product contains higher levels of dithiodicarbamate than natural latex rubber and consequently may cause sensitization by skin contact and give rise to Type IV Allergies. Persons who know they are susceptible to dithiocarbamates should avoid contact with polyisoprene. N-Nitrosamine by-products of dithiocarbamate are below BSEN12868 limits.
The substance is not classified, according to the GHS regulation.*
- **Label elements**
- **GHS label elements** *Void*
- **Hazard pictograms** *Void*
- **Signal word** *Void*
- **Hazard statements** *Void*

AU

(Contd. on page 2)

Safety Data Sheet

according to WHS Regulations

Printing date 28.07.2023

Revision: 28.07.2023

Trade name: **DermaDam™ Synthetic**

(Contd. of page 1)

3 Composition and Information on Ingredients

- **Chemical characterisation:** Mixtures
- **Description:** Mixture of substances listed below with nonhazardous additions.
9003-31-0 Polyisoprene Rubber

4 First Aid Measures

- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** If skin irritation continues, consult a doctor.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** None
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire Fighting Measures

- **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Protective equipment:** No special measures required.

6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and Storage

- **Handling:**
- **Precautions for safe handling:**
Avoid heat, sparks, and flames
Keep away from heat and direct sunlight.
No special measures required.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**
Store in dry conditions.
See product labelling.
Store in a cool place.

(Contd. on page 3)

Safety Data Sheet

according to WHS Regulations

Printing date 28.07.2023

Revision: 28.07.2023

Trade name: **DermaDam™ Synthetic**

(Contd. of page 2)

· **Specific end use(s)** Professional Dental Dam

8 Exposure controls and personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.
- **Ingredients with limit values that require monitoring at the workplace:** Not required.
- **Additional information:** The lists valid during the making were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:** Not required.
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
- **Penetration time of glove material**
The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Not required.
- **Body protection:** Protective work clothing

9 Physical and Chemical Properties

- **General Information**
- **Appearance:**
- **Form:** Solid
- **Colour:** Light blue
- **Odour:** Product specific
- **Odour threshold:** Not determined.
- **pH-value:** Not applicable.
- **Change in condition**
- **Melting point/freezing point:** Undetermined.
- **Initial boiling point and boiling range:** Undetermined.
- **Flash point:** Not applicable.
- **Flammability (solid, gas):** Product is not flammable.
- **Decomposition temperature:** Not determined.
- **Auto-ignition temperature:** Not determined.
- **Explosive properties:** Product does not present an explosion hazard.
- **Explosion limits:**
- **Lower:** Not determined.
- **Upper:** Not determined.
- **Vapour pressure:** Not applicable.
- **Density:** Not determined.
- **Relative density** Not determined.
- **Vapour density** Not applicable.
- **Evaporation rate** Not applicable.

(Contd. on page 4)

Safety Data Sheet

according to WHS Regulations

Printing date 28.07.2023

Revision: 28.07.2023

Trade name: **DermaDam™ Synthetic**

(Contd. of page 3)

- | | |
|---|---|
| <ul style="list-style-type: none"> · Solubility in / Miscibility with · water: · Partition coefficient: n-octanol/water: · Viscosity: · Dynamic: · Kinematic: | <p><i>Insoluble.</i></p> <p><i>Not determined.</i></p> <p><i>Not applicable.</i></p> <p><i>Not applicable.</i></p> |
| <ul style="list-style-type: none"> · Other information | <p><i>Partly soluble in organic solvents.</i></p> <p><i>Tensile Strength: 12.0 MPa</i></p> <p><i>Modulus at 500% Extension: 1.3 MPa</i></p> <p><i>Elongation at Break: 1100%</i></p> <p><i>Tear Strength: 24 N/mm</i></p> |

10 Stability and Reactivity

- **Reactivity** *No further relevant information available.*
- **Thermal decomposition / conditions to be avoided:** *No decomposition if used according to specifications.*
- **Possibility of hazardous reactions:** *No dangerous reactions known.*
- **Conditions to avoid:**
 - Ignition sources*
 - Sparks.*
 - Heat*
 - Flames*
- **Incompatible materials:**
 - Copper & Copper containing alloys*
 - Fats & Oils*
 - Solvents*
- **Hazardous decomposition products:** *No dangerous decomposition products known.*
- **Additional information:**

The behavioural properties of polyisoprene rubber with respect to its performances when exposed to chemicals are comparable to that of natural rubber.

11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** *Based on available data, the classification criteria are not met.*
- **Skin corrosion/irritation** *Based on available data, the classification criteria are not met.*
- **Serious eye damage/irritation** *Based on available data, the classification criteria are not met.*
- **Respiratory or skin sensitisation** *Based on available data, the classification criteria are not met.*
- **Germ cell mutagenicity** *Based on available data, the classification criteria are not met.*
- **Carcinogenicity** *Based on available data, the classification criteria are not met.*
- **Reproductive toxicity** *Based on available data, the classification criteria are not met.*
- **STOT-single exposure** *Based on available data, the classification criteria are not met.*
- **STOT-repeated exposure** *Based on available data, the classification criteria are not met.*
- **Aspiration hazard** *Based on available data, the classification criteria are not met.*

12 Ecological Information

- **Toxicity**
- **Aquatic toxicity:** *No further relevant information available.*
- **Persistence and degradability** *No further relevant information available.*

(Contd. on page 5)

Safety Data Sheet

according to WHS Regulations

Printing date 28.07.2023

Revision: 28.07.2023

Trade name: **DermaDam™ Synthetic**

(Contd. of page 4)

- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** *Disposal must be made according to official regulations.*

14 Transport information

· UN-Number	
· ADG, ADN, IMDG, IATA	not regulated
· UN proper shipping name	
· ADG, ADN, IMDG, IATA	not regulated
· Transport hazard class(es)	
· ADG, ADN, IMDG, IATA	
· Class	not regulated
· Packing group	
· ADG, IMDG, IATA	not regulated
· Environmental hazards:	Not applicable.
· Special precautions for user	Not Applicable
· Transport in bulk according to Annex II of Marpol and the IBC Code	Not applicable.
· UN "Model Regulation":	not regulated

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· Australian Inventory of Industrial Chemicals

Substance is listed.

· Standard for the Uniform Scheduling of Medicines and Poisons

Substance is not listed.

(Contd. on page 6)

Safety Data Sheet

according to WHS Regulations

Printing date 28.07.2023

Revision: 28.07.2023

Trade name: **DermaDam™ Synthetic**

(Contd. of page 5)

· **Australia: Priority Existing Chemicals**

Substance is not listed.

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** Substance is not listed.

· **Chemical safety assessment:** A chemical safety assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:** Environmental, Health, and Safety

· **Contact:** Customer Service

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

· *** Data compared to the previous version altered.**

AU