



Safety Data Sheet according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

*

1 Identification

- **Other means of identification**
- **Trade name:** UltraSeal XT™ Plus
- **Article number:**
SDS 25-001.18R01, 13560, 56303, 10952, 73401, 13526, 14907, 35202, 35203, 14869, 726, 727, 7661, 734, 1289, 565
- **Relevant identified uses of the substance or mixture and uses advised against**
Professional Dental Pit and Fissure Sealant
- **Application of the substance / the mixture** Professional Dental Pit and Fissure Sealant
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Ultradent Products, Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
(800) 552-5512

Ultradent Australia Pty Ltd.
Level 22/2 Market Street
Sydney NSW 2000
Australia
Email: info.anz@ultradent.com
Toll Free: 1-800-290929
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887

2 Hazard(s) Identification

- **Classification of the substance or mixture**



GHS07

Skin sensitisation – Category 1 H317 May cause an allergic skin reaction.

- **Additional information:** The full text of the classification and labelling elements mentioned is given in Section 16.

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard-determining components of labelling:**
Triethylene Glycol Dimethacrylate (≥10-<30 %)
Diurethane Dimethacrylate (≥1-<10 %)
- **Hazard statements**
H317 May cause an allergic skin reaction.
- **Precautionary statements**
P101 If medical advice is needed, have product container or label at hand.

(Contd. on page 2)

Safety Data Sheet

according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

Trade name: UltraSeal XT™ Plus

(Contd. of page 1)

- P102 Keep out of reach of children.
 P103 Read label before use.
 P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
 P280 Wear protective gloves.
 P362+P364 Take off contaminated clothing and wash it before reuse.
 P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
 P321 Specific treatment (see on this label).
 P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

3 Composition and Information on Ingredients

- **Chemical characterisation: Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

- **Dangerous components:**

109-16-0	Triethylene Glycol Dimethacrylate	≥10-<30%
	⚠ Skin sensitisation – Category 1, H317	
72869-86-4	Diurethane Dimethacrylate	≥1-<10%
	⚠ Skin sensitisation – Category 1, H317	
868-77-9	2-Hydroxyethyl Methacrylate	≥0.1-<1%
	⚠ Skin corrosion/irritation – Category 2, H315; Eye damage/irritation – Category 2A, H319; Skin sensitisation – Category 1, H317	
	Trade Secret	≥0.1-<1%
	⚠ Acute toxicity - oral – Category 4, H302; Acute toxicity - dermal – Category 4, H312; Skin corrosion/irritation – Category 2, H315; Eye damage/irritation – Category 2A, H319; Skin sensitisation – Category 1, H317; Flammable liquids – Category 4, H227	

- **Additional information:** For the wording of the listed hazard phrases refer to section 16.

4 First Aid Measures

- **General information:** Immediately remove any clothing soiled by the product.
- **After inhalation:**
 Supply fresh air and to be sure call for a doctor.
 In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If swallowed in large quantities seek medical attention.
- **Information for doctor:**
 · **Most important symptoms and effects, both acute and delayed** No further relevant information available.
 · **Indication of any immediate medical attention and special treatment needed**
 No further relevant information available.

5 Fire Fighting Measures

- **Suitable extinguishing agents:**
 Foam, dry chemical, carbon dioxide
 Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.

(Contd. on page 3)

Safety Data Sheet

according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

Trade name: UltraSeal XT™ Plus

(Contd. of page 2)

· **Protective equipment:** No special measures required.

6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures** Not required
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Dispose contaminated material as waste according to section 13.
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and Storage

- **Handling:**
- **Precautions for safe handling:**
Avoid contact with eyes, skin, and clothing. Use suitable protective equipment.
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** See product labeling.
- **Information about storage in one common storage facility:** Not required
- **Further information about storage conditions:** See product labelling.
- **Specific end use(s)** Professional Dental Pit and Fissure Sealant

8 Exposure controls and personal protection

- **Appropriate engineering controls** No further data; see section 7.
- **Ingredients with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists valid during the making were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
- **Respiratory protection:**
In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.
- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

(Contd. on page 4)

Safety Data Sheet

according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

Trade name: UltraSeal XT™ Plus

(Contd. of page 3)

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

- **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

- **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **Eye protection:** Goggles recommended during refilling

- **Body protection:** Protective work clothing

9 Physical and Chemical Properties

- **General Information**

- **Appearance:**

- **Form:** Medium-High Viscosity Liquid

- **Colour:** Shade dependent

- **Odour:** Acrylic

- **Odour threshold:** Not determined

- **pH-value:** Not applicable (non-aqueous)

- **Change in condition**

- **Melting point/freezing point:** Undetermined

- **Initial boiling point and boiling range:** Undetermined

- **Flash point:** >100 °C

- **Flammability** Not applicable

- **Decomposition temperature:** Not determined

- **Ignition temperature:** Product is not selfigniting.

- **Explosive properties:** Product does not present an explosion hazard.

- **Explosion limits:**

- **Lower:** Not determined

- **Upper:** Not determined

- **Vapour pressure at 20 °C:** 0 hPa

- **Density at 20 °C:** 1.7-1.8 g/cm³

- **Relative density** Not determined

- **Vapour density** Not determined

- **Evaporation rate** Not determined

- **Solubility in / Miscibility with**

- **water:** Partly miscible.

- **Partition coefficient: n-octanol/water:** Not determined

- **Viscosity:**

- **Dynamic:** Not determined

- **Kinematic:** Not determined

- **Other information**

- **Particle characteristics** Not applicable

- **Physical state** Liquid

10 Stability and Reactivity

- **Reactivity** No further relevant information available.

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.

- **Possibility of hazardous reactions:** No dangerous reactions known.

(Contd. on page 5)

Safety Data Sheet

according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

Trade name: UltraSeal XT™ Plus

(Contd. of page 4)

- **Conditions to avoid:** Avoid light exposure
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· LD/LC50 values relevant for classification:

109-16-0 Triethylene Glycol Dimethacrylate

Oral	LD50	>5,000 mg/kg (rat)
	LC50 Fish	16.4 mg/l (Fish) (Toxicity to fish)
Dermal	LD50	>2,000 mg/kg (mouse)

72869-86-4 Diurethane Dimethacrylate

Oral	LD50	>5,000 mg/kg (rat)
------	------	--------------------

868-77-9 2-Hydroxyethyl Methacrylate

Oral	LD50	3,275 mg/kg (mouse)
		>5,000 mg/kg (rat)
	LC50 Fish	>100 mg/l (Fish)
Dermal	LD50	>3,000 mg/kg (rabbit)
	LC50(Daphnia magna)	24.1 mg/l (daphnia)

Trade Secret

Oral	LD50	1,550 mg/kg (rat)
	LC50 Fish	19 mg/l (Fish)
Dermal	LD50	2,000 mg/kg (rabbit)
Inhalative	LC50/4 h	96 mg/l (rat)

- **Primary irritant effect:**
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation** May cause an allergic skin reaction.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

12 Ecological Information

· Toxicity

· Aquatic toxicity:

109-16-0 Triethylene Glycol Dimethacrylate

EC50	>100 mg/kg (Algae)
Biodegradability	28 days (Aerobic) (Biodegradability testing)
Aqua toxicity	32 mg/l (daphnia) (No Observed Effect Concentration)

(Contd. on page 6)

AU

Safety Data Sheet

according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

Trade name: UltraSeal XT™ Plus

(Contd. of page 5)

72869-86-4 Diurethane Dimethacrylate

EC50 >0.6 mg/kg (Algae)

Biodegradability 28 days (Aerobic) (Biodegradability testing)

868-77-9 2-Hydroxyethyl Methacrylate

EC50 345 mg/kg (Algae)

Trade Secret

EC50 42 mg/kg (Algae)

- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
Must not be disposed together with household garbage. Do not allow product to reach sewage system.
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

- **UN-Number**
- **ADG, ADN, IMDG, IATA** not regulated
- **UN proper shipping name**
- **ADG, ADN, IMDG, IATA** not regulated
- **Transport hazard class(es)**
- **ADG, ADN, IMDG, IATA**
- **Class** not regulated
- **Packing group**
- **ADG, IMDG, IATA** not regulated
- **Environmental hazards:** Not applicable
- **Special precautions for user** Not Applicable
- **Transport in bulk according to Annex II of Marpol and the IBC Code** Not applicable

(Contd. on page 7)

Safety Data Sheet

according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

Trade name: UltraSeal XT™ Plus

(Contd. of page 6)

· **UN "Model Regulation":** not regulated

*

15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **Australian Inventory of Industrial Chemicals**

65997-17-3 Inorganic Glass

109-16-0 Triethylene Glycol Dimethacrylate

72869-86-4 Diurethane Dimethacrylate

Trade Secret

868-77-9 2-Hydroxyethyl Methacrylate

64-17-5 Ethyl Alcohol

2530-85-0 Silane

10163-15-2 Sodium Monofluorophosphate

10373-78-1 Camphorquinone

Organophosphine Oxide

1332-37-2 Red Iron Oxide

51274-00-1 Yellow Iron Oxide

Chromium Iron Oxide

· **Standard for the Uniform Scheduling of Medicines and Poisons**

868-77-9 2-Hydroxyethyl Methacrylate

S5

· **Australia: Priority Existing Chemicals**

None of the ingredients is listed.

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** None of the ingredients is listed.

· **Chemical safety assessment:**

Device is biocompatible when used as directed by dental professionals per ISO 10993-1

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases from Section 3**

H227 Combustible liquid.

H302 Harmful if swallowed.

H312 Harmful in contact with skin.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

· **Department issuing SDS:** Regulatory Compliance

· **Contact:** Customer Service

(Contd. on page 8)

AU

Safety Data Sheet

according to WHS Regulations

Date of issue: 29.07.2026

Revision: 29.07.2026

Trade name: UltraSeal XT™ Plus

(Contd. of page 7)

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

Flammable liquids – Category 4: Flammable liquids – Category 4

Acute toxicity - oral – Category 4: Acute toxicity – Category 4

Skin corrosion/irritation – Category 2: Skin corrosion/irritation – Category 2

Eye damage/irritation – Category 2A: Serious eye damage/eye irritation – Category 2A

Skin sensitisation – Category 1: Skin sensitisation – Category 1

· *** Data compared to the previous version altered.**

AU