



Safety Data Sheet according to WHS Regulations

Date of issue: 30.04.2026



Revision: 30.04.2026

1 Identification

- **Other means of identification**
- **Trade name:** Consepsis™ Scrub
- **Article number:** SDS 403-001.04R01, 10969, 14493, 689
- **Relevant identified uses of the substance or mixture and uses advised against**
Professional Dental Chlorhexidine Gluconate Slurry
- **Application of the substance / the mixture** Professional Dental Chlorhexidine Gluconate Slurry
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Ultradent Products, Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
(800) 552-5512

Ultradent Australia Pty Ltd.
Level 22/2 Market Street
Sydney NSW 2000
Australia
Email: info.anz@ultradent.com
Toll Free: 1-800-290929
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887

2 Hazard(s) Identification

- **Classification of the substance or mixture**
-  GHS02 flame
Flammable liquids – Category 3 H226 Flammable liquid and vapour.
-  GHS07
Skin corrosion/irritation – Category 2 H315 Causes skin irritation.
Eye damage/irritation – Category 2A H319 Causes serious eye irritation.
- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** GHS02, GHS07
- **Signal word** Warning
- **Hazard statements**
H226 Flammable liquid and vapour.
H315 Causes skin irritation.
H319 Causes serious eye irritation.

(Contd. on page 2)

Safety Data Sheet

according to WHS Regulations

Date of issue: 30.04.2026

Revision: 30.04.2026

Trade name: **Consepsis™ Scrub**

(Contd. of page 1)

· Precautionary statements

- P101 *If medical advice is needed, have product container or label at hand.*
- P102 *Keep out of reach of children.*
- P103 *Read label before use.*
- P210 *Keep away from heat/sparks/open flames/hot surfaces. No smoking.*
- P241 *Use explosion-proof electrical/ventilating/lighting equipment.*
- P303+P361+P353 *IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.*
- P305+P351+P338 *IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.*
- P403+P235 *Store in a well-ventilated place. Keep cool.*
- P501 *Dispose of contents/container in accordance with local/regional/national/international regulations.*

3 Composition and Information on Ingredients

· Chemical characterisation: Mixtures

· **Description:** Mixture of substances listed below with nonhazardous additions.

· Dangerous components:

14808-60-7	Silica Glass ⚠ Acute toxicity - inhalation – Category 4, H332; Skin corrosion/irritation – Category 2, H315; Eye damage/irritation – Category 2A, H319	30-60%
64-17-5	Ethyl Alcohol ⚠ Flammable liquids – Category 2, H225; ⚠ Eye damage/irritation – Category 2A, H319	≥0-<10%
	Siloxane Polyalkyleneoxide Copolymer ⚠ Reproductive toxicity – Category 2, H361f	<3%
8006-90-4	Oils, Peppermint ⚠ Aspiration hazard – Category 1, H304; ⚠ Skin corrosion/irritation – Category 2, H315; Eye damage/irritation – Category 2A, H319; Skin sensitisation – Category 1, H317; Flammable liquids – Category 4, H227	≥0.1-<1%

· **Additional information:** For the wording of the listed hazard phrases refer to section 16.

4 First Aid Measures

- **General information:** Immediately remove any clothing soiled by the product.
- **After inhalation:** In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire Fighting Measures

· Suitable extinguishing agents:

CO₂ powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

(Contd. on page 3)

Safety Data Sheet

according to WHS Regulations

Date of issue: 30.04.2026

Revision: 30.04.2026

Trade name: **Consepsis™ Scrub**

(Contd. of page 2)

- **For safety reasons unsuitable extinguishing agents:** Water with full jet
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Protective equipment:** No special measures required.

6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures**
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Dispose contaminated material as waste according to section 13.
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and Storage

- **Handling:**
- **Precautions for safe handling:** No special precautions are necessary if used correctly.
- **Information about fire - and explosion protection:**
Keep ignition sources away - Do not smoke.
Protect against electrostatic charges.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required
- **Further information about storage conditions:**
See product labelling.
Keep container tightly sealed.
- **Specific end use(s)** Professional Dental Chlorhexidine Gluconate Slurry

8 Exposure controls and personal protection

- **Appropriate engineering controls** No further data; see section 7.
 - **Ingredients with limit values that require monitoring at the workplace:**
- | | |
|--------------------------------|------------------------------------------------------------|
| 14808-60-7 Silica Glass | |
| WES | Long-term value: 0.05 mg/m ³
respirable dust |
| 64-17-5 Ethyl Alcohol | |
| WES | Long-term value: 1880 mg/m ³ , 1000 ppm |
- **Additional information:** The lists valid during the making were used as basis.
 - **Personal protective equipment:**
 - **General protective and hygienic measures:**
Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.
 - **Respiratory protection:** Not required

(Contd. on page 4)

Safety Data Sheet

according to WHS Regulations

Date of issue: 30.04.2026

Revision: 30.04.2026

Trade name: **Consepsis™ Scrub**

(Contd. of page 3)

· **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**



Tightly sealed goggles

· **Body protection:** Protective work clothing

9 Physical and Chemical Properties

· **General Information**

· **Appearance:**

- **Form:** Paste
- **Colour:** Green
- **Odour:** Characteristic
- **Odour threshold:** Not determined
- **pH-value at 20 °C:** 5-7
- **Change in condition**
- **Melting point/freezing point:** Undetermined
- **Initial boiling point and boiling range:** Undetermined
- **Flash point:** 34 °C
- **Flammability** Flammable.
- **Auto-ignition temperature:** 363 °C
- **Decomposition temperature:** Not determined
- **Ignition temperature:** Product is not selfigniting.
- **Explosive properties:** Product is not explosive. However, formation of explosive air/vapour mixtures are possible.

· **Explosion limits:**

- **Lower:** 3.5 Vol %
- **Upper:** 15 Vol %
- **Vapour pressure:** Not determined
- **Density at 20 °C:** 1.57 g/cm³
- **Relative density** Not determined
- **Vapour density** Not determined
- **Evaporation rate** Not determined

(Contd. on page 5)

Safety Data Sheet

according to WHS Regulations

Date of issue: 30.04.2026

Revision: 30.04.2026

Trade name: **Consepsis™ Scrub**

(Contd. of page 4)

- | | |
|--------------------------------------------------|-----------------------------------|
| · Solubility in / Miscibility with | |
| · water: | Not miscible or difficult to mix. |
| · Partition coefficient: n-octanol/water: | Not determined |
| · Viscosity: | |
| · Dynamic: | Not determined |
| · Kinematic: | Not determined |
| · Other information | |
| · Particle characteristics | Not applicable |
| · Physical state | Liquid |

10 Stability and Reactivity

- **Reactivity** No further relevant information available.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:**
 - Ignition sources
 - Flames
 - Heat
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:**
 - Hydrogen chloride (HCl)
 - Oxides of Carbon
 - Nitrogen oxides
 - Ammonia

11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· LD/LC50 values relevant for classification:

ATE (Acute Toxicity Estimates)		
Inhalative	LC50/4 h	2.97 mg/l

64-17-5 Ethyl Alcohol

Oral	LD50	5,600 mg/kg (guinea pig)
		3,400 mg/kg (mouse)
		7,060 mg/kg (rat)
Inhalative	LC50 Fish	>10,000 mg/l (Fish)
	LC50/4 h	39 mg/l (mouse) 20,000 mg/l (rat)

8006-90-4 Oils, Peppermint

Oral	LD50	2,490 mg/kg (mouse)
		2,426 mg/kg (rat)

- **Primary irritant effect:**
- **Skin corrosion/irritation** Causes skin irritation.
- **Serious eye damage/irritation** Causes serious eye irritation.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.

(Contd. on page 6)

Safety Data Sheet

according to WHS Regulations

Date of issue: 30.04.2026

Revision: 30.04.2026

Trade name: **Consepsis™ Scrub**

(Contd. of page 5)

- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

12 Ecological Information

- **Toxicity**

- **Aquatic toxicity:**

64-17-5 Ethyl Alcohol	
------------------------------	--

Algae Toxicity	1,000 mg/l (Algae)
----------------	--------------------

- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
Disposal should be in accordance with applicable regional, national and local laws and regulations.
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** *Disposal must be made according to official regulations.*

14 Transport information

- | | |
|----------------------------------|---------------------------------------|
| · UN-Number | |
| · ADG, IMDG, IATA | UN1987 |
| · UN proper shipping name | |
| · ADG | 1987 ALCOHOLS, N.O.S. (Ethyl Alcohol) |
| · IMDG, IATA | ALCOHOLS, N.O.S. (Ethyl Alcohol) |

(Contd. on page 7)

AU

Safety Data Sheet

according to WHS Regulations

Date of issue: 30.04.2026

Revision: 30.04.2026

Trade name: **Consepsis™ Scrub**

(Contd. of page 6)

· **Transport hazard class(es)**· **ADG, IMDG, IATA**· **Class**

3 Flammable liquids.

· **Label**

3

· **Packing group**· **ADG, IMDG, IATA**

III

· **Environmental hazards:**

Not applicable

· **Special precautions for user**

Warning: Flammable liquids.

· **Hazard identification number (Kemler code):**

36

· **EMS Number:**

F-E,S-D

· **Stowage Category**

A

· **Transport in bulk according to Annex II of Marpol and the IBC Code**

Not applicable

· **Transport/Additional information:**· **ADG**· **Limited quantities (LQ)**

5L

· **Excepted quantities (EQ)**

Code: E1

Maximum net quantity per inner packaging: 30 ml

Maximum net quantity per outer packaging: 1000 ml

· **Transport category**

3

· **Tunnel restriction code**

D/E

· **IMDG**· **Limited quantities (LQ)**

5L

· **Excepted quantities (EQ)**

Code: E1

Maximum net quantity per inner packaging: 30 ml

Maximum net quantity per outer packaging: 1000 ml

· **UN "Model Regulation":**

UN 1987 ALCOHOLS, N.O.S. (ETHYL ALCOHOL), 3, III

15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

14808-60-7 Silica Glass

· **Australian Inventory of Industrial Chemicals**

14808-60-7 Silica Glass

7732-18-5 Water

25322-68-3 Polyethylene Glycol

64-17-5 Ethyl Alcohol

Trade Secret

18472-51-0 Chlorhexidine Gluconate

(Contd. on page 8)

Safety Data Sheet

according to WHS Regulations

Date of issue: 30.04.2026

Revision: 30.04.2026

Trade name: Consepsis™ Scrub

(Contd. of page 7)

22839-47-0	Aspartame
8006-90-4	Oils, Peppermint
5949-29-1	Citric Acid Monohydrate
89-78-1	Menthol

· Standard for the Uniform Scheduling of Medicines and Poisons

None of the ingredients is listed.

· Australia: Priority Existing Chemicals

None of the ingredients is listed.

· Directive 2012/18/EU
· Named dangerous substances - ANNEX I None of the ingredients is listed.

· Seveso category P5c FLAMMABLE LIQUIDS

· Qualifying quantity (tonnes) for the application of lower-tier requirements 5,000 t

· Qualifying quantity (tonnes) for the application of upper-tier requirements 50,000 t

· Chemical safety assessment:

Device is biocompatible when used as directed by dental professionals per ISO 10993-1

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Relevant phrases from Section 3

H225 Highly flammable liquid and vapour.

H227 Combustible liquid.

H304 May be fatal if swallowed and enters airways.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

H361f Suspected of damaging fertility.

· Department issuing SDS: Regulatory Compliance

· Contact: Customer Service

· Abbreviations and acronyms:

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

Flammable liquids – Category 2: Flammable liquids – Category 2

Flammable liquids – Category 3: Flammable liquids – Category 3

Flammable liquids – Category 4: Flammable liquids – Category 4

Acute toxicity - inhalation – Category 4: Acute toxicity – Category 4

Skin corrosion/irritation – Category 2: Skin corrosion/irritation – Category 2

Eye damage/irritation – Category 2A: Serious eye damage/eye irritation – Category 2A

Skin sensitisation – Category 1: Skin sensitisation – Category 1

Reproductive toxicity – Category 2: Reproductive toxicity – Category 2

Aspiration hazard – Category 1: Aspiration hazard – Category 1

· * Data compared to the previous version altered.