



Safety Data Sheet according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

1 Identification

- **Other means of identification**
- **Trade name:** VALO™ X Rechargeable Lithium-Ion Battery
- **Article number:** SDS 471-001.02R01, 5437
- **Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Battery
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Manufactured by GlobTek, Inc. for
Ultradent Products, Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
(800) 552-5512

Ultradent Australia Pty Ltd.
Level 22/2 Market Street
Sydney NSW 2000
Australia
Email: info.anz@ultradent.com
Toll Free: 1-800-290929
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887

2 Hazard(s) Identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonised System (GHS).
- **Additional information:**
Emergency Overview Safety Data Sheets (SDS) are a sub-requirement of the Occupational Safety and Health Administration (OSHA) Hazard Communication Standard, 29 CFR Subpart 1910.1200. This Hazard Communication Standard does not apply to various subcategories including anything defined by OSHA as an "article". OSHA has defined "article" as a manufactured item other than a fluid or particle; (i) which is formed to a specific shape or design during manufacture; (ii) which has end use function(s) dependent in whole or in part upon its shape or design during end use; and (iii) which under normal conditions of use does not release more than very small quantities, e.g. minute or trace amounts of hazardous chemical, and does not pose a physical hazard or health risk to employees. Because all of our batteries are defined as "articles", they are exempt from the requirements of the Hazard Communication Standard, hence an SDS is not required. However, this Safety Data Sheet (SDS) contains valuable information critical to the safe handling and proper use of this product. This SDS should be retained and available for employees and other users of this product.
- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void

(Contd. on page 2)

Safety Data Sheet

according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

· **Hazard statements** Void

(Contd. of page 1)

3 Composition and Information on Ingredients

- **Chemical characterisation:** Mixtures
- **Description:** Mixture of substances listed below with nonhazardous additions.

· **Dangerous components:**

7782-42-5	Graphite	10-<30%
21324-40-3	Lithium Hexafluorophosphate	10-<30%
	⚠ Acute toxicity - oral – Category 3, H301; Acute toxicity - dermal – Category 2, H310	
7440-50-8	Copper Foil	<10%
7440-02-0	nickel	≥0.1-<1%
	⚠ Carcinogenicity – Category 2, H351; Specific target organ toxicity (repeated exposure) – Category 1, H372; ⚠ Skin sensitisation – Category 1, H317	

· **Additional information:**

The chemicals listed in section 3 are contained in a sealed container. Risk of exposure only occurs if battery is mechanically, thermally, or electrically abused.

For the wording of the listed hazard phrases refer to section 16.

4 First Aid Measures

· **General information:**

The following information applies if the battery is mechanically, thermally or electrically abused.

· **After inhalation:**

If symptoms are experienced, remove source of contamination or move victim to fresh air. Get medical attention.

· **After skin contact:**

Flush affected area with lukewarm water for at least 30 minutes. If skin irritation persists, call a physician.

· **After eye contact:**

Immediately flush eyes with water for 30 minutes while lifting the upper and lower lids. Get medical attention.

· **After swallowing:**

Do not induce vomiting. Call a physician or Poison Control Center. National battery ingestion hotline: 202-625-333

· **Information for doctor:** Treat symptomatically.

· **Most important symptoms and effects, both acute and delayed**

Chemicals may cause burns to skin, eyes, gastrointestinal tract and mucous membranes. Contact with skin may cause chronic eczema or nickel itch. Electrolyte is extremely corrosive to eye tissue and may cause permanent blindness. If swallowed it may cause choking, nausea, persistent vomiting, diarrhea, abdominal pain, dizziness, faintness, unconsciousness, and possible liver and kidney injury.

· **Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

5 Fire Fighting Measures

· **Suitable extinguishing agents:**

Water spray

Foam

Dry powder

· **Special hazards arising from the substance or mixture**

Cells may rupture when exposed to excessive heat. This could result in the release of flammable or corrosive materials. Cells or batteries may flame or leak potentially hazardous organic vapors if exposed to excessive heat or fire. Damaged or opened cells or batteries can result in rapid heating and the release of flammable vapors.

(Contd. on page 3)

Safety Data Sheet

according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

(Contd. of page 2)

Vapors may be heavier than air and may travel along the ground or be moved by ventilation to an ignition source and flash back. LiPF salt contained in the electrolyte released hydrogen fluoride in contact with water.

· **Advice for firefighters:**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear. During water application, caution is advised as burning pieces of flammables particles may be ejected from the fire.

· **Protective equipment:** No special measures required.

· **Additional information Hazardous Combustion Products:** Carbon oxides. Hydrogen fluoride. Phosphorous oxides.

6 Accidental Release Measures

· **Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment as required. Wear protective gloves. Ventilate affected area.

· **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.

· **Methods and material for containment and cleaning up:**

Prevent further leakage or spillage if safe to do so.

Damaged batteries that are NOT hot or burning should be placed in a sealed plastic bag or plastic-lined metal container. If cells rupture and a thermal event follows: using a shovel or broom, cover battery or spilled substances with dry sand or vermiculite, place in approved container (after cooling if necessary).

· **Reference to other sections**

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and Storage

· **Handling:**

· **Precautions for safe handling:**

Do not expose battery or cell to extreme temperatures or fire. Do not disassemble, crush or puncture battery. Avoid mechanical or electrical abuse. Do not short circuit. Use only approved charges and charging procedures. Do not disassemble a battery any safety device.

· **Information about fire - and explosion protection:** No special measures required.

· **Storage:**

· **Requirements to be met by storerooms and receptacles:**

Batteries should be separated from other materials and stored in a non-combustible, well-ventilated, sprinkler-protected structure with sufficient clearance between walls and battery sacks. Do not place batteries near heating equipment; do not expose to direct sunlight for extended periods. Do not store batteries above 60 degrees Celsius or below -32 degrees Celsius. Store batteries in a cool (below 21 degrees Celsius or 70 degrees Fahrenheit), dry area that is subject to little temperature change. Elevated temperatures can result in reduced battery service life. Battery exposure to temperatures in excess of 130 degrees Celsius will result in the battery venting flammable liquid and gases. Do not store batteries in a manner that allows terminals to short circuit.

· **Information about storage in one common storage facility:** Not required.

· **Further information about storage conditions:**

Packaging Materials: If packing materials are not available, place masking tape on positive and negative ends of the cells.

Incompatible Materials: If leaked, forbidden to contact with strong oxidizers, mineral acids, strong alkalis, halogenated hydrocarbons. Water with internal contents of battery.

· **Specific end use(s)** No further relevant information available.

AU

(Contd. on page 4)

Safety Data Sheet

according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

(Contd. of page 3)

8 Exposure controls and personal protection

· **Appropriate engineering controls** No further data; see section 7.

· **Ingredients with limit values that require monitoring at the workplace:**

7782-42-5 Graphite

WES Long-term value: 3 mg/m³

7440-50-8 Copper Foil

WES Long-term value: 1 * 0.2** mg/m³
*dust & mists (as Cu) **fume

7440-02-0 nickel

WES Long-term value: 1 mg/m³
Metal: Sen

· **Additional information:** The lists valid during the making were used as basis.

· **Personal protective equipment:**

· **General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

· **Respiratory protection:**

Not required under normal conditions.

In case of battery venting or rupture, use a self contained full face respiratory mask. Refer to 29 CFR 1910.134 for respiratory protection requirements.

· **Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:** Safety glasses

· **Body protection:**

Protective work clothing

In case of battery rupture or leakage, wear rubber apron and Viton rubber gloves, Protective clothing. Refer to 29 CFR 1910.138 for appropriate skin and body protection.

9 Physical and Chemical Properties

· **General Information**

· **Appearance:**

· **Form:**

Geometric, solid object

· **Colour:**

According to product specification

· **Odour:**

Not determined

· **Odour threshold:**

Not determined.

· **pH-value:**

Not applicable.

· **Change in condition**

· **Melting point/freezing point:**

Undetermined.

· **Initial boiling point and boiling range:** Undetermined.

(Contd. on page 5)

Safety Data Sheet

according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

(Contd. of page 4)

· Flash point:	Not applicable.
· Flammability	Not determined.
· Decomposition temperature:	Not determined.
· Ignition temperature:	Product is not selfigniting.
· Explosive properties:	Product does not present an explosion hazard.
· Explosion limits:	
· Lower:	Not determined.
· Upper:	Not determined.
· Vapour pressure:	Not applicable.
· Density:	Not determined.
· Relative density	Not determined.
· Vapour density	Not applicable.
· Evaporation rate	Not applicable.
· Solubility in / Miscibility with	
· water:	Insoluble.
· Partition coefficient: n-octanol/water:	Not determined.
· Viscosity:	
· Dynamic:	Not applicable.
· Kinematic:	Not applicable.
· Other information	
· Particle characteristics	Not determined.
· Physical state	Solid

10 Stability and Reactivity

- **Reactivity** Not reactive under normal conditions of use.
- **Chemical stability** Stable under recommended conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** Reacts with water.
- **Conditions to avoid:**
Heat
Mechanical and electrical abuse. Electrical shorting. Battery exposure to temperatures in excess of 130 degrees Celsius will result in the battery venting flammable liquid and gases.
- **Incompatible materials:**
If leaked, forbidden to contact with strong oxidizers, mineral acids, strong alkalis, halogenated hydrocarbons. Water with internal contents of battery.
- **Hazardous decomposition products:**
No dangerous decomposition products known.
Hydrogen fluoride
Phosphorous Oxides
Carbon oxides

11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· LD/LC50 values relevant for classification:

ATE (Acute Toxicity Estimates)

Oral	LD50	>500 mg/kg
Dermal	LD50	>250 mg/kg

(Contd. on page 6)

Safety Data Sheet

according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

(Contd. of page 5)

- **Primary irritant effect:**
- **Skin corrosion/irritation**
Irritant effect on the skin
Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation**
Corrosive to eyes and may cause severe damage including blindness
Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

12 Ecological Information

- **Toxicity**
- **Aquatic toxicity:** Very toxic to aquatic life with long lasting effects
- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water
Do not allow product to reach ground water, water course or sewage system.
Danger to drinking water if even small quantities leak into the ground.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

- | | |
|----------------------------------|----------------------------|
| · UN-Number | |
| · ADG, IMDG, IATA | UN3480 |
| · UN proper shipping name | |
| · ADG | 3480 LITHIUM ION BATTERIES |
| · IMDG, IATA | LITHIUM ION BATTERIES |

(Contd. on page 7)

Safety Data Sheet according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

(Contd. of page 6)

· **Transport hazard class(es)**· **ADG, IMDG**· **Class**

9 Miscellaneous dangerous substances and articles.

· **Label**

9A

· **IATA**· **Class**

9 Miscellaneous dangerous substances and articles.

· **Label**

9A

· **Packing group**· **ADG, IMDG, IATA**

not regulated

· **Environmental hazards:**· **Marine pollutant:**

Symbol (fish and tree)

· **Special marking (ADG):**

Symbol (fish and tree)

· **Special precautions for user**

Warning: Miscellaneous dangerous substances and articles.

· **Hazard identification number (Kemler code):**

-

· **EMS Number:**

F-A,S-I

· **Stowage Category**

A

· **Stowage Code**

SW19 For batteries transported in accordance with SP 376 or SP 377 Category C, unless transported on a short international voyage.

· **Transport in bulk according to Annex II of Marpol and the IBC Code**

Not applicable.

· **Transport/Additional information:**· **ADG**· **Limited quantities (LQ)**

0

· **Excepted quantities (EQ)**

Code: E0

Not permitted as Excepted Quantity

· **Transport category**

2

· **Tunnel restriction code**

E

· **IMDG**· **Limited quantities (LQ)**

0

· **Excepted quantities (EQ)**

Code: E0

Not permitted as Excepted Quantity

· **UN "Model Regulation":**

UN 3480 LITHIUM ION BATTERIES

AU

(Contd. on page 8)

Safety Data Sheet

according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

(Contd. of page 7)

15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

7440-02-0	nickel
-----------	--------

· **Australian Inventory of Industrial Chemicals**

12190-79-3	Lithium Cobaltate
------------	-------------------

7782-42-5	Graphite
-----------	----------

21324-40-3	Lithium Hexafluorophosphate
------------	-----------------------------

7440-50-8	Copper Foil
-----------	-------------

24937-79-9	Poly Vinylidene Fluoride (PVDF)
------------	---------------------------------

7440-02-0	nickel
-----------	--------

9002-86-2	polyvinyl chloride
-----------	--------------------

9002-88-4	Polyethylene low density
-----------	--------------------------

· **Standard for the Uniform Scheduling of Medicines and Poisons**

None of the ingredients is listed.

· **Australia: Priority Existing Chemicals**

None of the ingredients is listed.

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** None of the ingredients is listed.

· **Qualifying quantity (tonnes) for the application of lower-tier requirements** 100 t

· **Qualifying quantity (tonnes) for the application of upper-tier requirements** 200 t

· **Chemical safety assessment:** A chemical safety assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases from Section 3**

H301 Toxic if swallowed.

H310 Fatal in contact with skin.

H317 May cause an allergic skin reaction.

H351 Suspected of causing cancer.

H372 Causes damage to organs through prolonged or repeated exposure.

· **Department issuing SDS:** Environmental, Health, and Safety

· **Contact:** Customer Service

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

Acute toxicity - oral – Category 3: Acute toxicity – Category 3

Acute toxicity - dermal – Category 2: Acute toxicity – Category 2

(Contd. on page 9)

Safety Data Sheet
according to WHS Regulations

Date of issue: 17.12.2025

Revision: 04.02.2025

Trade name: VALO™ X Rechargeable Lithium-Ion Battery

(Contd. of page 8)

Skin sensitisation – Category 1: Skin sensitisation – Category 1

Carcinogenicity – Category 2: Carcinogenicity – Category 2

Specific target organ toxicity (repeated exposure) – Category 1: Specific target organ toxicity (repeated exposure) – Category 1

· *** Data compared to the previous version altered.**

-AU-