



Safety Data Sheet according to WHS Regulations

Date of issue: 11.08.2026

Revision: 11.08.2026

*

1 Identification

- **Other means of identification**
- **Trade name:** Opalescence Go™ 15% Mint - Finished Good
- **Article number:**
SDS 478-1000.02R01, 1005955, 4633-US, 4638-AU, 4638-U, 4638-US, 4648-U, 4648-US, 6227-US, 6238-US, 6288-US, 1007781, 4894-US, S4633-AU, S4633-US
- **Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the mixture** Professional Dental Teeth Whitening Gel
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Ultradent Products, Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
(800) 552-5512

Ultradent Australia Pty Ltd.
Level 22/2 Market Street
Sydney NSW 2000
Australia
Email: info.anz@ultradent.com
Toll Free: 1-800-290929
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887

2 Hazard(s) Identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonised System (GHS).
- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void

*

3 Composition and Information on Ingredients

- **Chemical characterisation:** Mixtures
- **Description:** Mixture of substances listed below with nonhazardous additions.

(Contd. on page 2)

Safety Data Sheet

according to WHS Regulations

Date of issue: 11.08.2026

Revision: 11.08.2026

Trade name: Opalescence Go™ 15% Mint - Finished Good

(Contd. of page 1)

· Dangerous components:

56-81-5	Glycerin	<10%
	⚠ Eye damage/irritation – Category 2A, H319	
7722-84-1	Hydrogen Peroxide	3.69%
	⚠ Oxidising liquids - Category 1, H271; ⚠ Skin corrosion/irritation – Category 1A, H314;	
	⚠ Acute toxicity - oral – Category 4, H302; Acute toxicity - inhalation – Category 4, H332	

· **Additional information:** For the wording of the listed hazard phrases refer to section 16.

4 First Aid Measures

- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire Fighting Measures

- **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Protective equipment:** No special measures required.

6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures** Not required
- **Environmental precautions:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and Storage

- **Handling:**
- **Precautions for safe handling:**
See product labeling.
No special measures required.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required

(Contd. on page 3)

Safety Data Sheet

according to WHS Regulations

Date of issue: 11.08.2026

Revision: 11.08.2026

Trade name: Opalescence Go™ 15% Mint - Finished Good

(Contd. of page 2)

- **Further information about storage conditions:** None.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls and personal protection

- **Appropriate engineering controls** No further data; see section 7.

- **Ingredients with limit values that require monitoring at the workplace:**

56-81-5 Glycerin

TWA	Short-term value: 10 mg/m ³
WEL	Long-term value: 10 mg/m ³
WES	Long-term value: 10 mg/m ³ inhalable dust

7722-84-1 Hydrogen Peroxide

WEL	Long-term value: 0.7 mg/m ³ , 0.5 ppm
WES	Long-term value: 1.4 mg/m ³ , 1 ppm

- **Additional information:** The lists valid during the making were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:** Not required
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.
- **Penetration time of glove material**
The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Goggles recommended during refilling
- **Body protection:** Protective work clothing

9 Physical and Chemical Properties

- **General Information**
- **Appearance:**
- **Form:** Gel
- **Colour:** White Opaque
- **Odour:** Mint
- **Odour threshold:** Not determined
- **pH-value at 20 °C:** 4.8-6.8
- **Change in condition**
- **Melting point/freezing point:** Undetermined
- **Initial boiling point and boiling range:** Undetermined
- **Flash point:** Not applicable
- **Flammability** Not applicable

(Contd. on page 4)

AU

Safety Data Sheet

according to WHS Regulations

Date of issue: 11.08.2026

Revision: 11.08.2026

Trade name: Opalescence Go™ 15% Mint - Finished Good

(Contd. of page 3)

· Decomposition temperature:	Not determined
· Ignition temperature:	Product is not selfigniting.
· Explosive properties:	Product does not present an explosion hazard.
· Explosion limits:	
· Lower:	Not determined
· Upper:	Not determined
· Vapour pressure:	Not determined
· Density at 20 °C:	1.3 g/cm ³
· Relative density	Not determined
· Vapour density	Not determined
· Evaporation rate	Not determined
· Solubility in / Miscibility with	
· water:	Partly miscible.
· Partition coefficient: n-octanol/water:	Not determined
· Viscosity:	
· Dynamic:	Not determined
· Kinematic:	Not determined
· Other information	
· Particle characteristics	Not applicable
· Physical state	Liquid

10 Stability and Reactivity

- **Reactivity** No further relevant information available.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

LD/LC50 values relevant for classification:

ATE (Acute Toxicity Estimates)

Oral	LD50	26,269 mg/kg
Inhalative	LC50/4 h	578 mg/l

56-81-5 Glycerin

Oral	LD50	7,750 mg/kg (guinea pig) 4,100 mg/kg (mouse) 5,570 mg/kg (rat) 27,000 mg/kg (rabbit)
	LC50 Fish	>5,000 mg/l (Fish)
Dermal	LD50	>21,900 mg/kg (rat) 10,000 mg/kg (rabbit)

(Contd. on page 5)

Safety Data Sheet

according to WHS Regulations

Date of issue: 11.08.2026

Revision: 11.08.2026

Trade name: Opalescence Go™ 15% Mint - Finished Good

(Contd. of page 4)

7722-84-1 Hydrogen Peroxide

Oral	LC50 Fish	16.4 mg/l (Fish)
------	-----------	------------------

- **Primary irritant effect:**
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

12 Ecological Information· **Toxicity**· **Aquatic toxicity:****56-81-5 Glycerin**

EC50	>10,000 mg/kg (Bacteria)
------	--------------------------

7722-84-1 Hydrogen Peroxide

EC50	1.38 mg/l (Algae)
	2.4 mg/l (daphnia)

- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
 Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
 Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
 Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

AU

(Contd. on page 6)

Safety Data Sheet

according to WHS Regulations

Date of issue: 11.08.2026

Revision: 11.08.2026

Trade name: Opalescence Go™ 15% Mint - Finished Good

(Contd. of page 5)

14 Transport information

· UN-Number	
· ADG, IMDG, IATA	not regulated
· UN proper shipping name	
· ADG, ADN, IMDG, IATA	not regulated
· Transport hazard class(es)	
· ADG, ADN, IMDG, IATA	
· Class	not regulated
· Packing group	
· ADG, IMDG, IATA	not regulated
· Environmental hazards:	Not applicable
· Special precautions for user	Not Applicable
· Transport in bulk according to Annex II of Marpol and the IBC Code	Not applicable
· UN "Model Regulation":	not regulated

15 Regulatory information

· Safety, health and environmental regulations/legislation specific for the substance or mixture		
· NIOSH-Ca (National Institute for Occupational Safety and Health)		
None of the ingredients is listed.		
· Australian Inventory of Industrial Chemicals		
56-81-5	Glycerin	
7732-18-5	Water	
7722-84-1	Hydrogen Peroxide	
9003-39-8	Polyvinylpyrrolidone	
1310-58-3	Potassium Hydroxide	
	Trade Secret	
7758-11-4	Dipotassium Phosphate	
87-99-0	Xylitol	
7757-79-1	Potassium Nitrate	
	Citric Acid, Anhydrous	
8006-90-4	Oils, Peppermint	
7681-49-4	Sodium Fluoride	
56038-13-2	Sucralose	
6381-92-6	Edetate Disodium	
68585-47-7	Sodium Lauryl Sulfate	
79-10-7	acrylic acid	
110-82-7	cyclohexane	
· Standard for the Uniform Scheduling of Medicines and Poisons		
7722-84-1	Hydrogen Peroxide	S5, S6, S10
1310-58-3	Potassium Hydroxide	S5, S6, S10

(Contd. on page 7)

Safety Data Sheet

according to WHS Regulations

Date of issue: 11.08.2026

Revision: 11.08.2026

Trade name: Opalescence Go™ 15% Mint - Finished Good

(Contd. of page 6)

· Australia: Priority Existing Chemicals
None of the ingredients is listed.
· Directive 2012/18/EU
· Named dangerous substances - ANNEX I *None of the ingredients is listed.*
· Chemical safety assessment: *A chemical safety assessment has not been carried out.*

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Relevant phrases from Section 3
H271 May cause fire or explosion; strong oxidizer.
H302 Harmful if swallowed.
H314 Causes severe skin burns and eye damage.
H319 Causes serious eye irritation.
H332 Harmful if inhaled.
· Department issuing SDS: *Regulatory Compliance*
· Contact: *Customer Service*
· Abbreviations and acronyms:
ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
IATA: International Air Transport Association
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
LC50: Lethal concentration, 50 percent
LD50: Lethal dose, 50 percent
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
Oxidising liquids - Category 1: Oxidizing liquids – Category 1
Acute toxicity - oral – Category 4: Acute toxicity – Category 4
Skin corrosion/irritation – Category 1A: Skin corrosion/irritation – Category 1A
Eye damage/irritation – Category 2A: Serious eye damage/eye irritation – Category 2A
· * Data compared to the previous version altered.

AU