



Safety Data Sheet according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

*

1 Identification

- **Other means of identification**
- **Trade name:** Opalescence™ Custom 35% Birthday Cake
- **Article number:** SDS 534-001.01R01, 7047-US, 1014097, 7179-US
- **Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Ultradent Products, Inc.
505 W. Ultradent Drive (10200 S)
South Jordan, UT 84095-3942
USA
onlineordersupport@ultradent.com
(800) 552-5512

Ultradent Australia Pty Ltd.
Level 22/2 Market Street
Sydney NSW 2000
Australia
Email: info.anz@ultradent.com
Toll Free: 1-800-290929
- **Further information obtainable from:** Customer Service
- **Emergency telephone number:**
CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300
(INTERNATIONAL) : +(703) 527-3887

2 Hazard(s) Identification

- **Classification of the substance or mixture**



- Skin corrosion/irritation – Category 2 H315 Causes skin irritation.
Eye damage/irritation – Category 2A H319 Causes serious eye irritation.
· **Additional information:** The full text of the classification and labelling elements mentioned is given in Section 16.

- **Label elements**
- **GHS label elements** The product is classified and labelled according to the Globally Harmonised System (GHS).
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard statements**
H315 Causes skin irritation.
H319 Causes serious eye irritation.
- **Precautionary statements**
P101 If medical advice is needed, have product container or label at hand.
P102 Keep out of reach of children.
P103 Read label before use.
P264 Wash thoroughly after handling.

(Contd. on page 2)

AU

Safety Data Sheet

according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 1)

- P280 Wear protective gloves / eye protection / face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P332+P313 If skin irritation occurs: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before reuse.
P337+P313 If eye irritation persists: Get medical advice/attention.

3 Composition and Information on Ingredients

- **Chemical characterisation: Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

- **Dangerous components:**

56-81-5	Glycerin	≥10-<30%
	⚠ Eye damage/irritation – Category 2A, H319	
124-43-6	Carbamide Peroxide	10-<30%
	⚠ Oxidising solids - Category 3, H272; ⚠ Skin corrosion/irritation – Category 1B, H314; Eye damage/irritation – Category 1, H318	
7722-84-1	Hydrogen Peroxide	≥5-<8%
	⚠ Oxidising liquids - Category 1, H271; ⚠ Skin corrosion/irritation – Category 1A, H314; ⚠ Acute toxicity - oral – Category 4, H302; Acute toxicity - inhalation – Category 4, H332	
1310-73-2	Sodium Hydroxide	≥1-<3%
	⚠ Acute toxicity - oral – Category 3, H301; ⚠ Skin corrosion/irritation – Category 1A, H314; Eye damage/irritation – Category 1, H318; ⚠ Acute toxicity - dermal – Category 4, H312	

- **Additional information:** For the wording of the listed hazard phrases refer to section 16.

4 First Aid Measures

- **General information:** Immediately remove any clothing soiled by the product.
- **After inhalation:** In case of unconsciousness place patient stably in side position for transportation.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:**
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire Fighting Measures

- **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Protective equipment:** No special measures required.

6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures** Not required

(Contd. on page 3)

Safety Data Sheet

according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 2)

- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**
 See Section 7 for information on safe handling.
 See Section 8 for information on personal protection equipment.
 See Section 13 for disposal information.

7 Handling and Storage

- **Handling:**
- **Precautions for safe handling:** See product labeling.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required
- **Further information about storage conditions:** Keep container tightly sealed.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls and personal protection

- **Appropriate engineering controls** No further data; see section 7.
- **Ingredients with limit values that require monitoring at the workplace:**

56-81-5 Glycerin

TWA	Short-term value: 10 mg/m ³
WEL	Long-term value: 10 mg/m ³
WES	Long-term value: 10 mg/m ³ inhalable dust

7722-84-1 Hydrogen Peroxide

WEL	Long-term value: 0.7 mg/m ³ , 0.5 ppm
WES	Long-term value: 1.4 mg/m ³ , 1 ppm

1310-73-2 Sodium Hydroxide

WEL	Peak limitation: 2 mg/m ³
WES	Peak limitation: 2 mg/m ³

- **Additional information:** The lists valid during the making were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**
 Keep away from foodstuffs, beverages and feed.
 Immediately remove all soiled and contaminated clothing.
 Wash hands before breaks and at the end of work.
 Avoid contact with the eyes and skin.
- **Respiratory protection:** Not required
- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
 Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

(Contd. on page 4)

Safety Data Sheet

according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 3)

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· **Material of gloves**

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:**



Tightly sealed goggles

· **Body protection:** Protective work clothing

9 Physical and Chemical Properties

· **General Information**

· **Appearance:**

· **Form:**

Gel

· **Colour:**

Clear

· **Odour:**

Birthday Cake

· **Odour threshold:**

Not determined

· **pH-value at 20 °C:**

5.4-7.4

· **Change in condition**

· **Melting point/freezing point:**

Undetermined

· **Initial boiling point and boiling range:**

Undetermined

· **Flash point:**

Not applicable

· **Flammability**

Not determined

· **Decomposition temperature:**

Not determined

· **Ignition temperature:**

Product is not selfigniting.

· **Explosive properties:**

Product does not present an explosion hazard.

· **Explosion limits:**

· **Lower:**

1.5 Vol %

· **Upper:**

15 Vol %

· **Vapour pressure:**

Not applicable

· **Vapour pressure at 50 °C:**

280 hPa

· **Density:**

Not determined

· **Relative density**

Not determined

· **Vapour density**

Not applicable

· **Evaporation rate**

Not applicable.

· **Solubility in / Miscibility with**

· **water:**

Soluble.

· **Partition coefficient: n-octanol/water:**

Not determined

· **Viscosity:**

· **Dynamic:**

Not applicable

· **Kinematic:**

Not applicable

· **Other information**

· **Particle characteristics**

Not determined

(Contd. on page 5)

AU

Safety Data Sheet

according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 4)

· **Physical state** Solid

10 Stability and Reactivity

- **Reactivity** No further relevant information available.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· LD/LC50 values relevant for classification:

ATE (Acute Toxicity Estimates)

Oral	LD50	3,002-5,120 mg/kg
Dermal	LD50	46,552 mg/kg
Inhalative	LC50/4 h	200 mg/l

56-81-5 Glycerin

Oral	LD50	7,750 mg/kg (guinea pig) 4,100 mg/kg (mouse) 5,570 mg/kg (rat) 27,000 mg/kg (rabbit)
Dermal	LC50 Fish	>5,000 mg/l (Fish)
	LD50	>21,900 mg/kg (rat) 10,000 mg/kg (rabbit)
Inhalative	LC50/4 h	>0.1425 mg/l (rat)

124-43-6 Carbamide Peroxide

Oral	LD50	>2,000 mg/kg (rat)
------	------	--------------------

7722-84-1 Hydrogen Peroxide

Oral	LC50 Fish	16.4 mg/l (Fish)
------	-----------	------------------

1310-73-2 Sodium Hydroxide

Oral	LD50	130-340 mg/kg (rat)
	LC50 Fish	160 mg/l (Fish)
Dermal	LD50	1,350 mg/kg (rabbit)
	Absolute lethal concentration	180 ppm (Fish)

- **Primary irritant effect:**
- **Skin corrosion/irritation** Causes skin irritation.
- **Serious eye damage/irritation** Causes serious eye irritation.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.

(Contd. on page 6)

Safety Data Sheet

according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 5)

- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

12 Ecological Information

- **Toxicity**

- **Aquatic toxicity:**

- 56-81-5 Glycerin

- EC50 >10,000 mg/kg (Bacteria)

- 7722-84-1 Hydrogen Peroxide

- EC50 1.38 mg/l (Algae)

- 2.4 mg/l (daphnia)

- 1310-73-2 Sodium Hydroxide

- EC50 40.38 mg/kg (Water Flea)

- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
 - Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water
 - Do not allow product to reach ground water, water course or sewage system.
 - Danger to drinking water if even small quantities leak into the ground.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
 - Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

14 Transport information

- **UN-Number**
- **ADG, IMDG, IATA** not regulated
- **UN proper shipping name**
- **ADG, IMDG, IATA** not regulated

(Contd. on page 7)

Safety Data Sheet

according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 6)

· Transport hazard class(es)	
· ADG, ADN, IMDG, IATA	
· Class	not regulated
· Packing group	
· ADG, IMDG, IATA	not regulated
· Environmental hazards:	Not applicable
· Special precautions for user	Not Applicable
· Transport in bulk according to Annex II of Marpol and the IBC Code	Not applicable
· UN "Model Regulation":	not regulated

15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **Australian Inventory of Industrial Chemicals**

All ingredients are listed.

· **Standard for the Uniform Scheduling of Medicines and Poisons**

124-43-6 Carbamide Peroxide

S5, S6, S10

7722-84-1 Hydrogen Peroxide

S5, S6, S10

· **Australia: Priority Existing Chemicals**

None of the ingredients is listed.

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** None of the ingredients is listed.

· **Chemical safety assessment:** A chemical safety assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases from Section 3**

H271 May cause fire or explosion; strong oxidizer.

H272 May intensify fire; oxidizer.

H301 Toxic if swallowed.

H302 Harmful if swallowed.

H312 Harmful in contact with skin.

H314 Causes severe skin burns and eye damage.

H318 Causes serious eye damage.

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

· **Department issuing SDS:** Regulatory Compliance

· **Contact:** Customer Service

(Contd. on page 8)

Safety Data Sheet

according to WHS Regulations

Date of issue: 14.09.2026

Revision: 14.09.2026

Trade name: Opalescence™ Custom 35% Birthday Cake

(Contd. of page 7)

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

Oxidising liquids - Category 1: Oxidizing liquids – Category 1

Oxidising solids - Category 3: Oxidizing solids – Category 3

Acute toxicity - oral – Category 3: Acute toxicity – Category 3

Acute toxicity - oral – Category 4: Acute toxicity – Category 4

Skin corrosion/irritation – Category 1A: Skin corrosion/irritation – Category 1A

Skin corrosion/irritation – Category 1B: Skin corrosion/irritation – Category 1B

Skin corrosion/irritation – Category 2: Skin corrosion/irritation – Category 2

Eye damage/irritation – Category 1: Serious eye damage/eye irritation – Category 1

Eye damage/irritation – Category 2A: Serious eye damage/eye irritation – Category 2A

· *** Data compared to the previous version altered.**

AU